

Manfred Sauer Nephrostomy Codes

Description	Inlet length / tap type	Qty	Order Code	PIP Code	Select code required. 'X' (Healthcare Professional (HCP) to select)
500 ml sterile drainage bags only, 20cm inlet, luer lock connector	20cm / Lever tap	10 bags only	NS721.1720S	305-1547	
500 ml sterile drainage bags only, 30 cm inlet, luer lock connector	30 cm / Lever tap	10 bags only	NS721.1730S	305-1539	
500 ml sterile drainage bags only, 20cm inlet, luer lock connector	20 cm / Twist tap	10 bags only	NS721.3720S	399-1585	
500 ml sterile drainage bags only, 30 cm inlet, luer lock connector	30 cm / Twist tap	10 bags only	NS721.3730S	399-1593	
Replacement belt - adult	-	1 belt	NSBelt.02	305-1513	
Replacement belt - small	-	1 belt	NSBelt.01	305-1521	
Night Bag - 2L capacity drainable (for NephSys)	120 cm / Lever tap	10 bags only	NB2	320-3569	
Night Bag - 3L capacity drainable (for NephSys)	130 cm / Lever tap	10 bags only	NB3D	400-4610	
Day Bag to Night Bag Connectors - Latex-free	-	10 pcs	55.39-10	384-7514	
P.Grip Tabs x 10, 1 x hook-and-loop strip	-	1 pc	CS.TABS	304-6570	
Tegaderm dressing	15cm x20cm	1 pc	1628	003-8703	
Drainguard catheter securement device code		1 pc	8150		
Respond adhesive remover spray		50ml	OPN50	390-8845	

*PIP codes help pharmacy identify products.

To Register, Call

0800 18 19 01

Email Us

sales@rapidcare.co.uk



Lever Tap



Twist Tap